



An overview of coherent X-ray microscopy techniques

Ana Diaz

Paul Scherrer Institute (Switzerland)

Seminar

ALBA, 7th April 2025

The cSAXS team



Mads A.
Carlsen



Ana
Diaz




Xavier
Donath



Kazu
Hirosawa




Mirko
Holler



Andreas
Menzel



Clemens
Schmid



Postdoc
position

Alumni

I. Rodriguez-Fernandez	K. Wakonig	J. Han
T. Aidukas	D. Karpov	O. Bunk
N. Phillips	M. Odstrcil	H. Deyhle
C. Appel	V. Lütz-Bueno	T. Ikonen
M. Guizar-Sicairos	E. Tsai	C. Kewish
J. Ihli	M. Liebi	F. Pfeiffer
Z. Gao	I. Rajkovic	P. Thibault
M. Verezhak	R. Jacob	M. Dierolf
	J. da Silva	

Collaborations within PSI

J. Raabe (CPS)
C. David (CPS)
E. Föjd (detectors)
M. Guizar-Sicairos (CPS & EPFL)
G. Aeppli (CPS)
M. Liebi (CPS & EPFL)
K. Wakonig (CSD)
C. Appel (CSD)

S. Shahmoradian (CLS)
T. Ishikawa (CLS)
E. Müller (CLS)
A. Wanner (CLS)
B. Rösner (CPS)
H.-C. Stadler (CSD)
L. Debenjak (CAS)

External collaborators (PI's)

M.A.G. Aranda (Uni. Malaga, Spain)
C. Bosch (Francis Crick Inst, UK)
A. Schaefer (Francis Crick Inst, UK)
A. Pacureanu (ESRF, France)

Introduction

- Motivation to phase contrast hard X-ray microscopy
- Basic concepts: X-ray wave propagation, coherence...

Phase-contrast hard X-ray microscopy



- Zernike phase-contrast microscopy
- Holo-tomography
- Coherent diffraction imaging
- Ptychography

Outlook

- High-brilliance synchrotron sources


Conventional X-ray microscopy

full-field microscopy:




M. Uchida *et al.*,
Yeast **28**, 227 (2011)

scanning transmission microscopy (STXM):




B. Watts *et al.*,
Materials Today
15, 148 (2012)



- Typically with Fresnel zone plates in the soft X-ray regime (e.g. in the water window)
- Resolution can reach 25-30 nm, but depth of focus limited to about 1 μm


Imaging bulk samples with hierarchical structures

Computer chip



<https://www.livescience.com/52207-faster-3d-computer-chip.html>

Nervous tissue



Copyright © 2007 Pearson Education, Inc., publishing as Benjamin Cummings


Thickness from 10 to 100 μm
Resolution from 10 to 100 nm

Energy > 2 keV


Challenges:

- Low absorption
- Fabrication aberration-free, high-resolution lenses

Hard X-ray transmission microscopy with absorption contrast



S. Spence *et al.*, Nanotechnol. **32**, 442003 (2021)



V. De Andrade *et al.*, Adv. Mater. **33**, 2008653 (2021)

Absorption and phase contrast






Image from “Phase-contrast X-ray imaging” in Wikipedia:
https://en.wikipedia.org/wiki/Phase-contrast_X-ray_imaging

For hard X-rays: $\delta \gg \beta$

Full-field microscopy with Zernike phase contrast

M. Stampanoni *et al.*, Phys. Rev. B **81**, 140105(R) (2010)



- Holes matching custom illumination
- The transmitted beam goes through the holes
- The beam which is refracted by the sample goes through the phase-shift mask
- Both parts contribute to a phase-contrast image on the detector
- Low coherence requirement. Problems if the beam has a high degree of coherence

Wave propagation reveals phase contrast





Figure adapted from P. Willmott, "An Introduction to Synchrotron Radiation", Wiley (2011) p. 314

Transversal coherence

J. F. Van der Veen & F. Pfeiffer
J. Phys.: Condens. Matter **16**, 5003 (2004)




Transversal coherence length: $\xi_t = \frac{\lambda R}{w}$

Maximum distance between two slits such that they produce constructive interference when illuminated by an extended source size

Transversal coherence

J. F. Van der Veen & F. Pfeiffer
J. Phys.: Condens. Matter **16**, 5003 (2004)



Transversal coherence length:


$$\xi_t = \frac{\lambda R}{w}$$

$$\lambda = 1 \text{ \AA}$$

$$R = 50 \text{ m}$$

$$w_v = 20 \text{ \mu m} \rightarrow \xi_v = 250 \text{ \mu m}$$


$$w_h = 200 \text{ \mu m} \rightarrow \xi_h = 25 \text{ \mu m}$$



- Interference between refracted wave and incoming beam: aka **in-line holography**
- Intensity fringes build up in the near field, especially at the edges of the sample

(*) See e.g. P. C. Diemoz *et al.*, Opt. Express **20**, 2789 (2012)


Propagation-based phase contrast




P. Cloetens *et al.*,
PNAS **103**, 14626 (2006)

Polystyrene foam
 $\lambda = 0.69 \text{ \AA}$


$d = 0.01 \text{ m}$




$d = 0.12 \text{ m}$



$d = 0.22 \text{ m}$



$d = 0.91 \text{ m}$




P. Cloetens *et al.*, J. Phys D: Appl. Phys **32**, A145 (1999)


100 μm

Propagation-based phase contrast: holotomography


Quantitative phase image
obtained with a
reconstruction algorithm



Tomographic reconstruction




Slice through tomographic
reconstruction



P. Cloetens *et al.*, Appl. Phys. Lett. **75** 1912 (1999)


P. Cloetens *et al.*, J. Phys D: Appl. Phys **32**, A145 (1999)

Propagation-based imaging with magnification



- Focused beam creates divergent illumination onto the specimen
- Resolution limited to focus size
- Similar algorithms as in propagation-based phase-contrast imaging
- Equivalent propagation distance D and magnification M :

$$D = \frac{D_1 D_2}{D_1 + D_2} \quad M = \frac{D_1 + D_2}{D_1}$$



A. T. Kuan et al.,
Nat. Neurosci. **23**, 1637-1643 (2020)


M. Langer *et al.*, PLOS ONE **7**, e35691 (2012)

R. Mokso *et al.*, Appl. Phys. Lett. **90**, 144104 (2007)


Also known as **holo-nanotomography**

Lens- or focusing-based hard X-ray microscopy methods


Full-field microscopy with Zernike phase-contrast
M. Stampanoni *et al.*, Phys. Rev. B **81**, 140105(R) (2010)



Nano-holotomography
M. Langer *et al.*, PLOS ONE **7**, e35691 (2012)




Wave propagation reveals phase contrast




Coherent diffraction imaging (CDI)

support constraint




Sample consisting of
100 nm Au nanodots

Iterative phase retrieval algorithm



J.R. Fienup Appl. Opt. **21**, 2758 (1956)



oversampling



- Resolution not limited by optics


J. Miao et al., Nature **400**, 342 (1999)

Coherent diffraction imaging (CDI) on single cells






H. Jiang *et al.*,
PNAS **107**, 11234 (2010)

- Metal-coated polymer micro-sphere




E. T. B. Skjønsvell *et al.*,
J. Opt. Soc. Am. A **35**, 7-17 (2018)

Bragg coherent diffraction imaging (Bragg CDI)




3D displacement field in nano-crystals




M. Pfeifer *et al.*, Nature **442** (2006) 63

CDI: practical limitations




J.R. Fienup,
Appl. Opt. **21**, 2758 (1956)
J. Miao *et al.*,
Nature **400**, 342 (1999)

- Nyquist sampling \Rightarrow confined sample
- Typical beamline setups
 - Sample-detector ~ 5 m
 - Pixel size ~ 50 μm
 - Wavelength ~ 1 -2 \AA \Rightarrow Sample size $<$ a few micron
- Convergence




These limitations can be overcome using a modulator or a structured illumination

Coherent diffraction patterns from overlapping illuminated areas




- Absorption and phase contrast
- Resolution not limited by a lens!
- In practice limited by mechanical stability and thermal drifts

Iterative phase retrieval algorithms to reconstruct complex-valued transmissivity




H. M. L. Faulkner & J. M. Rodenburg,
Phys. Rev. Lett. **93**, 023903 (2004)

A typical X-ray ptychography setup




J. Vila-Comamala *et al.*, Opt. Express **19**, 21333 (2011)



coherent flux:
 5×10^8 photons/s
@ 6.2 keV


Simultaneous probe reconstruction with ptychography

Phase image





8 nm resolution
2×2 μm area
4 min


Illumination beam



Ptychographic X-ray computed tomography (PXCT)



Mouse bone specimen



M. Dierolf *et al.*, Nature **467** (2010) 436

Voxel size: 65 nm
Resolution: 120 nm
Dose: 2MGy


OMNY: tOMography Nano crYo stage


M. Holler and J. Raabe

- Laser interferometry for relative positioning of sample and illumination optics
- Aimed 3D resolution: 10 nm
- Cryo stage in ultra-high vacuum
- First test setup in air at room temperature, still in user operation

M. Holler *et al.*, Rev. Sci. Instrum. **83**, 073703 (2012)

M. Holler *et al.*, Rev. Sci. Instrum. **89**, 043706 (2018)






- Robust algorithms for online processing
- **Sample needs to be surrounded by air on both sides at all angles**
- Tomographic reconstructions are provided to the user during the experiment

M. Guizar-Sicairos *et al.*,
Opt. Express **19**, 21345 (2011)
M. Odstrčil *et al.*,
Opt. Express **27**, 36637 (2019)


Detection of synapses in stained, resin-embedded tissue



3D electron density map



38 nm resolution
absorbed dose: 2.5×10^9 Gy



Adrian
Wanner
PSI



Tomas
Aidukas
PSI



Carles
Bosch
Francis Crick
Institute

C. Bosch *et al.*,
bioRxiv (2024) doi: 10.1101/2023.11.16.567403
(in review)

- Novel sample preparation with radiation-hard resin
- Non-rigid tomographic reconstruction for sample deforming during acquisition




PXCT on 7 nm-node chip

4 nm 3D resolution
Burst ptychography
Record in hard X-ray microscopy

T. Aidukas *et al.*, Nature **632**, 81 (2024)




Tomas Aidukas




CDI vs ptychography

coherent diffraction imaging



H. Jiang *et al.*, PNAS **107**, 11234 (2010)


ptychography



M. Guizar-Sicairos *et al.*, Opt. Express **19**, 21345 (2011)



- 3D reconstruction from 3D reciprocal space: very robust to positioning accuracy
- Limited to isolated samples of a few microns
- Requires cutting-edge instrumentation for high-resolution scanning + rotation
- Applicable to large samples

Near-field ptychography



Nano-holotomography

M. Langer *et al.*, PLOS ONE **7**, e35691 (2012)





Near-field ptychography



M. Stockmar *et al.*, Sci. Rep. **3**, 1927 (2013)

Comparing PXCT and nano-holotomography


ptycho-tomography
@ 6.2 keV
108 nm, 3.3e7 Gy



holo-tomography
@ 17 keV
121 nm, 8.6e7 Gy



Heavy-metal-stained brain tissue




Alexandra
Pacureanu
ESRF


- Similar spatial resolution ~ 100 nm
- 2.5x less dose with ptycho-tomo
- 16x faster with holo-tomography (in resels/s)

Collaboration with
A. Pacureanu (ESRF, France)
C. Bosch & A. Schaefer (Francis Crick Inst., UK)
Manuscript in review



The upgrade: SLS 2.0 and cSAXS 2.0



- ~ 30x more coherent flux with SLS 2.0
- ~ 3x more flux with new U17 undulator
- ~ 10x more flux with more efficient optics for PXCT
- ~ 40x more flux with broadband option
- Better beam stability



Benefits of upgraded source



coherent flux:
 5×10^{10} photons/s
@ 6.2 keV

J. Vila-Comamala *et al.*, Opt. Express **19**, 21333 (2011)




Radiation damage
Detector saturation

Proposals to overcome the challenge of high photon flux

- I. Scan faster ☐ careful, same flux density (photons/s/μm²)
- II. Spread flux on sample and on detector ☐ no high-resolution scanning probe possible

Far/near-field ptycho
holotomography
high energy




large
single photon counting
detector

CoDi 😊

Conclusions – high-resolution hard X-ray microscopy

Transmission X-ray microscopy (TXM)


Very fast measurements
Inherently dose inefficient
Resolution limited by lens



S. Spence *et al.*,
Nanotechnol. **32**, 442003
(2021)

X-ray holo-nanotomography (XNH)


Fast measurements
Phase contrast
Resolution limited by focus



A. T. Kuan *et al.*, Nat.
Neurosci. **23**, 1637-1643
(2020)

Coherent diffraction imaging (CDI)


Confined, small sample
Convergence issues
Ongoing developments



H. Jiang *et al.*,
PNAS **107**, 11234 (2010)

X-ray ptychography

Extended samples
Good convergence
Challenging experiment



T. Aidukas *et al.*,
Nature **632**, 81 (2024)

Thank you for your attention
Questions?