## **REXS2025- International Conference on Resonant Elastic X-ray** Scattering

# **Beyond the Surface: 3D Probing of Antiferromagnets** A Journey Through Failures and Breakthroughs

## Valerio Scagnoli

Laboratory for Mesoscopic Systems, Department of Materials, ETH Zurich, 8093 Zurich, Switzerland PSI Center for Neutron and Muon Sciences, 5232 Villigen PSI, Switzerland

#### **ABSTRACT**

Three dimensional magnetic systems hold the promise to provide new functionality associated with greater degrees of freedom. Over the last years we have worked towards

developing methods to fabricate and characterize three dimensional magnetic structures. Specifically, we have combined X-ray magnetic imaging via circular dichroism (XMCD) with new iterative reconstruction algorithms to X-rav (ferro)magnetic tomography laminography in 3D volumes with sub 100 nm spatial [1-4]. Recent revival resolution of interest antiferromagnets have driven our recent efforts in developing an approach to image the antiferromagnetic order parameter in micron-size sample with nm spatial resolution. To this end, we have attempted to use coherent magnetic diffraction and antiferromagnetic tomographic imaging via X-ray linear dichroism. Over the course of several experimental campaigns, we have advanced these approaches and ultimately developed Xray Linear Dichroic Orientation Tomography (XL-DOT) a novel, quantitative, and non-invasive technique for threedimensional characterization of extended polycrystalline and non-crystalline materials at the intra- and intergranular The figure included here shows a [5-6].




Figure 1 Volume rendering of the c-axis orientation tomogram, obtained via XL-DOT, reveal the sample's microstructure. The coloured rods indicate the local alignment of the crystallographic c axis within each voxel, with colours representing the c-axis orientation projected onto the x-y plane.

reconstructed grain structure along with the local crystallographic c-axis alignment. The spectroscopic and non-destructive nature of XL-DOT makes it ideally suited for operando investigations, enabling simultaneous chemical and microstructural analysis of functional materials, including antiferromagnets.

### REFERENCES

- C. Donnelly et al., Nature 547, 328 (2017), <a href="https://doi.org/10.1038/nature23006">https://doi.org/10.1038/nature23006</a>.
  C. Donnelly et al., New J. Phys. 20, 083009 (2018), <a href="https://doi.org/10.1088/1367-2630/aad35a">https://doi.org/10.1088/1367-2630/aad35a</a>.
- 3. C. Donnelly et al., Nat. Phys. 17, 316 (2021), https://doi.org/10.1038/s41567-020-01057-3.
- 4. C. Donnelly et al., Nat. Nanotechnol. 15, 356 (2020), https://doi.org/10.1038/s41565-020-0649-x.
- 5. A. Apseros et al., Nature **636**, 354 (2024), <a href="https://doi.org/10.1038/s41586-024-08233-y">https://doi.org/10.1038/s41586-024-08233-y</a>.
- 6. A. Apseros et al., submitted to New J. Phys., https://arxiv.org/abs/2504.12978.